



Katarzyna Kamińska (1,2), Ewelina Nalejska (1,2), Magdalena Ratajska PhD (3), Bartosz Wasąg PhD(3), Tomasz Mierzwa MD PhD (4), Janusz Kowalewski MD PhD (1), Marzena Anna Lewandowska PhD MBA (1,2)

1. Department of Molecular Oncology and Genetics, Innovative Medical Forum, The F. Łukaszczyk Oncology Center, Bydgoszcz, Poland
2. Department of Thoracic Surgery And Tumors, The Ludwik Rydygier Collegium Medicum in Bydgoszcz, Nicolaus Copernicus University, Torun, Poland
3. Department of Biology and Medical Genetics, Medical University of Gdansk, Gdansk, Poland.
4. Department of Prevention and Health Promotion, The F. Łukaszczyk Oncology Center, Bydgoszcz, Poland



Contact: Marzena Anna Lewandowska Molecular Oncology and Genetics and Department, IFM, The Professor Franciszek Łukaszczyk Oncology Center. Romanowskiej 2, 85-796 Bydgoszcz, lewandowskam@co.bydgoszcz.pl

INTRODUCTION

Targeted therapy with PARP inhibitors inhibits growth of breast and ovarian cancer tumor cells harboring a BRCA1 or BRCA2 mutation. NGS analysis detects both somatic and germline mutations in tumors. Depending on the material (blood/tissue) and selected patient population (e.g. patients with platinum-sensitive ovarian cancer), mutations in the BRCA1/2 genes are found in 11-38% of ovarian and 7-9% breast cancer patients.

AIM

Evaluate the frequency of mutations in the BRCA1 and BRCA2 genes using NGS in a selected group of patients with breast or ovarian cancer.

MATERIALS AND METHODS

Patient eligibility criteria:

- Participated in the National Program for Combating Cancer (between 2006–2015);
- Negative results for the 3-5 most frequent BRCA1/2 mutations;
- Breast (n=45) or ovarian cancer (n=30).

1 Isolate DNA

2 DNA QC

3 Library prep

4 Library QC

5 Sequencing

6 Data analysis

DNA is isolated from FFPE tissue using QIAamp DNA FFPE Tissue Kit(Qiagen). Isolated DNA is qualified and quantified using Real-Time PCR by calculating the ratio of short to long DNA fragments (F-ratio) in a sample using the DNA Fragmentation Quantification Assay (EntroGen).

Libraries are made using BRCA MASTR™ Dx CE-IVD (Multiplicom) for ovarian cancer samples and BRCA Complete™ CE-IVD (EntroGen) for breast cancer samples. Libraries are qualified using the Library Quantification Kit for Illumina® (EntroGen).


Libraries are sequenced on the MiniSeq platform (Illumina). Bioinformatic analysis is performed under the Moncodaneum project (OncoDNA&Wallon Region, Belgium) using the OncoKDM tools developed by OncoDNA.

RESULTS

2. Sequencing the BRCA1 and BRCA2 coding regions of DNA isolated from tumor tissues found **7 additional mutations (23% more)** in ovarian cancer patients who were tested negative for BRCA1/2 mutations with the preventative screening program (2006-2015).

Table 2. List of mutations found in ovarian cancer patients (n=30). The first two samples marked with an asterisk (*) are positive controls with known mutations.

ID	Mutation	Classification	Consequence(s)
1.*	BRCA1 NM_007294.3:c.5266dupC (p.Gln1756Profs)	Pathogenic	Change in reading frame
2.*	BRCA1 NM_007294.3:c.181T>G (p.Cys61Gly)	Pathogenic	Amino acid change (missense mutation)
3.	BRCA1 NM_007294.3:c.850C>T (p.Gln284Ter)	Pathogenic	Amino acid change, STOP codon (nonsense mutation)
4.	BRCA1 NM_007294.3:c.406dupA (p.Arg136Lysfs)	Pathogenic	Change in reading frame
5.	BRCA1 NM_007294.3:c.4484+1G>A	Pathogenic	Amino acid change, splice donor site variant
6.	BRCA2 NM_000059.3:c.10095delCinsGAATTATAT CT (p.Ser3366Asnfs)	Likely pathogenic/ deleterious benign	Change in reading frame changes STOP codon (nonsense mutation), mRNA splicing (?) Literature: deleterious benign
7.	BRCA2 NM_000059.3:c.6267_6269delGCAinsC (p.Glu2089Aspfs)	Pathogenic	Change in reading frame
8.	BRCA2 NM_000059.3:c.2410G>T (p.Glu804Ter)	Pathogenic	Amino acid change, STOP codon (nonsense mutation)
9.	BRCA2 NM_000059.3:c.6267_6269delGCAinsC (p.Glu2089Aspfs)	Pathogenic	Change in reading frame




3. Mutations are absent in FFPE matched-normal tissue (marked in images above): BRCA1 c.5266dupC (1); BRCA1 c.181T>G (2) (capillary sequencer Applied Biosystems SeqStudio Genetic Analyzer); samples are negative for other mutations as well.




RESULTS

BREAST CANCER



OVARIAN CANCER



- The average age of diagnosis is 43.
 - 44% of patients were diagnosed at age 40 or earlier.
 - The youngest patient included in the study was diagnosed at age 25.
- The average age of diagnosis is 47.
 - 71% of patients were diagnosed at age 40 or earlier.
 - The youngest patient included in the study was diagnosed at age 25.

1. Sequencing the BRCA1 and BRCA2 coding regions of DNA isolated from tumor tissues found **4 additional mutations (about 9% more)** in breast cancer patients who were tested negative for BRCA1/2 mutations with the preventative screening program (2006-2015).

Table 1. List of mutations found in breast cancer patients (n=45).

ID	Mutation	Classification	Consequence(s)
1.	BRCA1 NM_007294.3:c.4689C>G (p.Tyr1563Ter)	Pathogenic	Amino acid change, STOP codon (nonsense mutation)
2.	BRCA1 NM_007294.3:c.5123C>A (p.Ala1708Glu)	Pathogenic	Amino acid change (missense mutation)
3.	BRCA1 NM_007294.3:c.1687C>T (p.Gln563Ter)	Pathogenic	Amino acid change, STOP codon (nonsense mutation)
4.	BRCA2 NM_000059.3:c.7758G>A (p.Trp2586Ter)	Pathogenic	Amino acid change, STOP codon (nonsense mutation)

CONCLUSIONS

- NGS makes it possible to detect both germline and somatic mutations in the BRCA1 and BRCA2 genes in breast and ovarian cancer tissues.
- Next generation sequencing increases the chance of detecting mutations in the BRCA1 and BRCA2 genes by 23% in patients with ovarian cancer and by 9% in patients with breast cancer compared to the national screening test for the 3-5 most common BRCA mutations in the Polish population.
- Sequencing a matched-normal tissue sample can determine whether a mutation is germline or somatic.
- Older FFPE material resulted in low uniformity of coverage due to DNA fragmentation. Therefore, there is a risk of false negative results in samples with a F-ratio <0.5.
- The detected mutations highlight the importance of bioinformatics tools used to analyze next-generation sequencing results.

Literature:

1. Ratajska M, Krygier M, Stukan M, et al. Mutational analysis of BRCA1/2 in a group of 134 consecutive ovarian cancer patients. Novel and recurrent BRCA1/2 alterations detected by next generation sequencing. Journal of Applied Genetics. 2015;56(2):193-198.
2. Koczkowska M, Zuk M, Gorczynski A, Ratajska M, Lewandowska M, Biernat W, Limon J, Wasag B. Detection of somatic BRCA1/2 mutations in ovarian cancer - next-generation sequencing analysis of 100 cases. Cancer Med. 2016 Jul;5(7):1640-6.
3. Dann RB, DeLoia JA, Timms KM, Zorn KK, Potter J, Flake DD 2nd, Lanchbury JS, Krivak TC. BRCA1/2 mutations and expression: response to platinum chemotherapy in patients with advanced stage epithelial ovarian cancer. Gynecol Oncol. 2012 Jun;125(3):677-82.